

## bisReadMapper\_v1 USER MANUAL

*This software only works on Linux operating systems (64 bits), and have been tested on:*

Ubuntu 10.04 LTS

Perl v5.10.1

SOAP 2.21release (Only works on 64 bits systems!)

Samtools 0.1.8

### Download required softwares:

1) SOAP - <http://soap.genomics.org.cn/download/soap2.21release.tar.gz>

After downloading, extract the files with:

```
tar -xvzf soap2.21release.tar.gz
```

2) Samtools - <http://sourceforge.net/projects/samtools/files/samtools/>

After downloading, extract the files with:

```
bzip2 -d samtools-0.1.8.tar.bz2
tar -xvf samtools-0.1.8.tar
```

Make sure that these dependencies are installed for samtools to compile:

```
sudo apt-get install zlib1g-dev libncurses5-dev
```

Then navigate into the folder and compile:

```
cd samtools-0.1.8
make
```

3) Check that Perl version 5 or later is installed

```
perl --version
```

### Download genome references:

1) full reference genome can be downloaded from UCSC Genome Browser. Note: chromosome names should be the standard names.

#### For hg18:

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg18/chromosomes/chr*.fa.gz
```

#### For hg19:

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg19/chromosomes/chr*.fa.gz
```

Move unwanted chromosome files from the directory before concatenating files to get standard chromosomes only.

```
i.e. mkdir nonStandardChromosomes
```

```
i.e. mv chr*_*.fa.gz > nonStandardChromosomes/
```

To get a single genome file, concatenate all chromosome files into one with less:

```
i.e. less chr*.fa.gz > HsGenome18.fa
```

For SNP calling, download the latest reference dbSNP file compatible with your specific reference genome:

#### For hg18:

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg18/database/snp130.txt.gz
```

#### For hg19:

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg19/database/snp132.txt.gz
```

To decompress dbSNP file:

```
i.e. gzip -d snp130.txt.gz
```

Required steps:

1) run genomePrep on single genome file (program will write to the current directory) . This program have two

optional inputs: context=[all, cg], and convert=[yes,no]. The default is to perform calls only on CpGs (context=cg, and convert=yes) and to perform in silico bisulfite conversion (generates \*.bis.fwd and \*.bis.rev files.) If convert=no is used, the converted files will not be generated.

i.e:     /path/to/genomePrep.pl genome.fasta  
Perform in silico bisulfite conversion and call only CpGs.

i.e:     /path/to/genomePrep.pl genome.fasta context=all convert=no  
Do not perform in silico bisulfite conversion and call all C (when you already have a converted genome.)

2) run 2bwt-builder (Go to <http://soap.genomics.org.cn/soapaligner.html> for more instructions) on the \*.bis.fwd and \*.bis.rev files.

i.e:     /path/to/soap/2bwt-builder genome.fasta.bis.fwd  
i.e:     /path/to/soap/2bwt-builder genome.fasta.bis.rev

3) run samtools faidx on the \*.bis.fwd and \*.bis.rev files

i.e:     /path/to/samtools/samtools faidx genome.fasta.bis.fwd  
i.e:     /path/to/samtools/samtools faidx genome.fasta.bis.rev

4) prepare params.txt file for each sample (See example file: sampleParams.txt)

<b>reads=</b>	provide full path to reads file(s), for pair-end, list two reads files separated by a comma (,) and no space in between
<b>name=</b>	(optional) provide the name of the sample to be used for labeling output files [default=Sample]
<b>allC=</b>	(optional) to call methylation on all C position (1) or only CG will be called (0) [default=0]
<b>length=</b>	(optional) provide the read length for calculating allowable mismatches (1 per 40 for single-end and 1 per 30 for paired-end) [default=80]
<b>refDir=</b>	provide path to reference genome index directory Note: there should only be one bisulfite reference per directory Required contents for refDir: *.bis.fwd, *.bis.rev, SOAP *.bis.fwd.index, SOAP *.bis.rev.index, *.chrSizes, *.cpositions.txt
<b>snp=</b>	(optional) provide full path to dbSNP file Note: No SNP will be called if this file is not provided
<b>soapDir=</b>	provide path to the SOAP software
<b>samtoolsDir=</b>	provide path to the Samtools software
<b>alignMode=</b>	specify whether to align reads as paired (P) or single-end (S) Note: must provide two reads files for setting P
<b>qualBase=</b>	Illumina base quality scale, Phred33 (33) or Phred64 (64)
<b>numCPU=</b>	(optional) specify the number of processor cores available for mapping [default:1]
<b>trim5=</b>	Number of base to trim from the 5' end of reads [default:0]
<b>trim3=</b>	Number of base to trim from the 3' end of reads [default:0]

4) run bisReadMapper.pl (program will output to the current directory)

i.e: `/path/to/bisReadMapper_v1.2/bisReadMapper.pl sampleParams.txt > sample.log`